

A trusted resource.

Since its release, the TotalCare® bed system has enabled critical care nurses from around the world to deliver patient care when it is most critical.

Throughout the years, proven performance has truly met continuous innovation as the TotalCare® bed continues to provide a consistent solution for both patients and their caregivers.



A trusted resource, at the most critical moment.

EASY PATIENT HANDLING

From lying in bed, to sitting upright, to standing up, it is important to have the right products in place to deliver care to your patients.

EASE OF USE AT EVERY POINT OF CARE

Patients in the ICU have unique medical needs that can make your routine more challenging. By using products that allow you to do more on your own, you have the flexibility to give patients the treatment they need – whenever they need it.

SAFE SKIN AT ANY LEVEL OF ACUITY

Due to the lack of relative mobility, patients in the ICU are susceptible to skin breakdown. As a result, it is important that their overall care prevents, treats and provides protection against pressure ulcers.



Progressive Mobility® Program

GUIDING PATIENTS ON THE PATH TO RECOVERY, ONE STEP AT A TIME.

Getting your patients moving and independent is one of the most challenging tasks you face. Working together with the TotalCare SpO₂RT® system and patient lifts, you can provide your patients with the support they need at the most critical moments to confidently progress from lying in the bed to sitting upright to standing up.



Head-of-bed and early Continuous Lateral Rotation Therapy





The tilt table position gradually acclimates patients to the upright position





The FullChair® position allows a single caregiver to more easily comply with an upright mobility protocol with the simple press of a button



↓↓↓ STAND

Chair egress and lifts permit safe, easy patient egress from the front





Patient lifts enable safe weight bearing exercise and/or walking support



Surfaces

Hill-Rom is a trusted partner at every level of acuity. From the non-powered NP100 foam surface to the powered air, weight-based pressure redistribution surface of the TotalCare SpO_2RT° 2 surface, the TotalCare SpO_2RT° 3 surface SpO_2RT° 4 surface SpO_2RT° 5 surface SpO_2RT°



TotalCare SpO₂RT® 2 Surface

- + Powered air, weight-based pressure redistribution surface
- + Prevention and treatment of skin complications associated with immobility
- + Comfort adjust options for patient satisfaction
- + Prevention of pressure ulcers in patients who are at moderate to high risk of skin breakdown
- + Treatment of pressure ulcers up to Stage III and IV
- + Therapeutic weight limit: 500 lbs



TotalCare® P500 Surface

- + Powered air, weight-based pressure redistribution surface
- + Prevention and treatment of skin complications associated with immobility
- + Prevention of pressure ulcers in patients who are at moderate to high risk of skin breakdown
- + Treatment of pressure ulcers up to Stage III and IV
- + Therapeutic weight limit: 500lbs



Accumax Quantum™ VPC/TC Surface

- + Non-powered Controlled Release™ Technology
- + Prevention and treatment of pressure ulcers up to and including Stage III and IV in patients with multiple turning surfaces
- + Weight limit: 500 lbs



TotalCare® NP100 Surface

- + Non-powered zoned foam sections allow for improved pressure redistribution points
- + Reduces the potential for pressure ulcers
- + Weight limit: 500 lbs

Ease of use at every point of care

Whether you're transporting patients or performing routine care, doing your job efficiently while also comfortably supporting your patient isn't easy.

Do more on your own with a range of features that help you turn, transport and reposition patients with a touch of a button—every step of the way.

- + IntelliDrive® Powered Transport allows a single caregiver to easily transport a patient in bed
- + Turn Assist allows you to turn a patient on your own, making it easier to perform linen changes or other routine tasks
- + A large, colorful touch screen offers control for functionality such as bed exit, integrated scale or head-of-bed alarms all at a touch of a button



Safe skin at any level of acuity

Keeping ICU patients' skin healthy throughout their recovery is one of your highest priorities. Achieve your goals and go beyond standard care with surfaces that address pressure, shear, friction and moisture—the four main causes of pressure ulcers.



PREVENT

FlexAfoot™ retractable foot mechanism makes it easy for you to adjust the overall length of the bed according to a patient's height—preventing patients from sliding down in bed and minimizing the need for repositioning.

PROTECT

Advanced Microclimate® technology helps reduce heat and moisture at the surface, keeping the patient's skin cool and dry.

TREAT

The pressure redistribution surface automatically adjusts to a patient's weight, body type and movement in both bed and FullChair® positioning to help treat and prevent pressure ulcers.

The TotalCare® Bed family: adapt to the needs of your individual patient.

SURFACE FEATURES

- + Advanced Microclimate® Technology helps reduce heat and moisture at the surface keeping the patient's skin cool and dry.
- + X-ray cassette sleeve built into the surface for easier cassette placement.
- + Automatic seat deflate during chair exit or manual activation for side of bed exit or bedpan placement.
- + Antimicrobial nano Ag+® technology integrated into surface kills up to 99.8% of stain and odor-causing bacteria.
- + Integrated air surface automatically adjusts according to patient's weight, body type, movement and bed position.

FRAME FEATURES

- + All siderails must be up to activate Turn Assist and CLRT.

 Lowering a siderail stops CLRT and alarms to alert caregiver.
- + One-step HandsFree® CPR max inflates the surface and flattens the frame, even from FullChair® position.
- + Siderails stow under bed frame for 'zero' transfer gap.
- + No siderail pinch-points protect fingers and critical lines from being pinched or pulled on when raising HOB.
- + No footend controls allow for easy removal of footboard for patient procedures, chair exit and cleaning.
- + Head-of-bed angle alarm and trending helps caregivers ensure protocol compliance.
- + Boost® feature max inflates the surface and flattens the frame with one button, with optional 7° Trendelenburg for easier repositioning.



CLRT OPTIONS

- + Full body continuous lateral rotation therapy (CLRT) maintains body alignment with optimal envelopment and pressure redistribution.
- + Percussion & Vibration (P&V) bladders directly under top of mattress for greater contact with patient.
- + Easy-to-use pre-set or customizable settings for CLRT and P&V to meet specific patient needs.

PATIENT MOBILITY

- + Cradle transition gently rocks frame back and slightly deflates surface during transition to FullChair® position for patient comfort and safety.
- + FlexAfoot™ feature built-in adjustable bed length expands and retracts bed length to fit patient height and for leg exercises.
- + Tilt Table position achieves up to 20° Reverse Trendelenburg for orthostatic and weight-bearing training.
- + One-button upright FullChair® position with head-of-bed at 65° and legs dependent at 70° with foot support.
- + Three options to get patients out of bed: side of bed, chair footend exit and chair side exit.
- + IntelliDrive® Powered Transport allows a single caregiver to



Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2012 Hill-Rom Services, Inc. ALL RIGHTS RESERVED. 177842 rev 1 08/03/2012 ENG – US

www.hill-rom.com

USA 800-445-3730 CANADA 800-267-2337 Enhancing outcomes for patients and their caregivers:

